



## SPECTROMETER

# *PRODUCT SAFETY CARD*

---

The Inphaze Spectrometer uses power from the power point the same way a normal PC Computer does.

## **IN CASE OF EMERGENCY:**

Push the Emergency Power Stop Button of your Workstation/Lab, if available – Please consult with your Lab Technician.

-- or --

Switch off the power point that supplies the Spectrometer's power, if available – Please consult with your Lab Technician.

-- or --

Switch off the Spectrometer's power switch. The Spectrometer's power switch is located at the rear, as shown below:

